

detect and identify

**NEW – NEW – NEW**

# CentroLIA<sub>pc</sub>

Microplate Luminometer



Berthold Technologies introduces

[CentroLIA<sub>pc</sub> LB 962](#)

## Microplate Luminometer



[Learn more about CentroLIA<sub>pc</sub>](#)

See us at



HALL 1, E 28

[Learn more about Luminescence Applications](#)

The CentroLIA<sub>pc</sub> LB 962 is a microplate luminometer unique in its design and performance:

- **crosstalk: 10<sup>-6</sup>**
- **detection limit: 0.001 mIU/L (TSH)**
- **user-friendly software**
- **immunoassay data reduction**



The [CentroLIA<sub>pc</sub>](#) is consequent development of the renowned CentroLIA. While maintaining the excellent performance the new CentroLIA<sub>pc</sub> has a smaller footprint and is more lightweight.

The instrument is controlled via an easy-to-use wizard driven PC software which includes immunassay data reduction.

The CentroLIA<sub>pc</sub> is an ideal instrument for diagnostic applications and clinical use as well as an attractive alternativ in research.

[Get in touch with Berthold Technologies](#)

[www.Berthold.com/Bio](http://www.Berthold.com/Bio)

To unsubscribe from this newsletter service please click here: [Bio.unsubscribe@berthold.com](mailto:Bio.unsubscribe@berthold.com)